

Eleganza 5



THE ART OF EARLY MOBILIZATION

MobiLift®

Simple to use sit-to-stand device.

MobiRails

Ergonomically shaped handles in siderails.

Lateral Tilt

Tilting improves the angle of the pelvis to assist in standing.

TURN LEFT, TURN RIGHT, DONE!

With frame-based lateral tilt of up to 15 degrees, the Eleganza 5 brings safety and comfort:



Eleganza 5 at a glance

- 1 Foot controls
Efficient hands-free control of the bed.
- 2 Electric Bed
Extension
- 3 Mattress platform
Easily removable and cleanable mattress platform.
- 4 Mechanical CPR
Electrical CPR activation on control panels.
- 5 5th castor
Improves maneuverability when transporting a patient.
- 6 Integrated patient controls
With illuminated keyboards.
- 7 Bed ends
Ergonomically shaped bed ends for easy transport.
- 8 iBrake
Automatic brake application.



USB Charging Port Standard



Optional iDRIVE Power available

TECHNICAL SPECIFICATIONS

Mattress dimensions	208 x 90 cm (81.9 x 35.4 in)	Max Backrest Angle	65°
Max. siderails height	45 cm (17.7 in)	Max. Bed Extension	22 cm / 8.7 in
Height adjustment	43.5 – 81.5 cm (17 x 32 in)		
Lateral tilt	± 15°	External dimensions (shortest position)	219x100 cm (86.2 x 39.4 in)
Safe working load	250 kg/551 Lbs	External dimensions (longest position)	241x100 cm (94.8 x 39.4 in)
Trendelenburg/Reverse Trendelenburg	± 14°	Frame Weight Std. Config.	210kg/463Lbs.